


CLARITY Tissue Clearing



Hydrogel Matrix
The first step in the tissue processing workflow is to place the fixed mouse, rat, or human tissue sample in a solution of hydrogel monomers and cross-linkers.



Functionalization
The monomers and cross-linkers diffuse into the tissue's cells and bind to biomolecules such as proteins and nucleic acids but not to light-scattering lipids.

Polymerization
The hydrogel is thermally treated and the monomers polymerize into stable mesh that locks proteins / biomolecules in place.

Lipid Clearing
A detergent removes lipids and other unbound molecules from the tissue. The proteins, nucleic acids and other bound biomolecules remain embedded with the hydrogel mesh.

Immunostaining
If desired, antibody-based immunostaining or labeling for many nucleic acids (RNA/DNA) can be used in a multiplex panel to highlight specific structures in the clarified sample.

Imaging
The tissue processing workflow moves from molecular labeling to imaging. The stained tissue is placed in a refractive index-matching mounting solution for imaging with a confocal or light sheet microscope or another imaging platform.



* Standard Resolution is dependent upon image resolution.
- improving resolution decreases the FOV
- better resolution equals smaller FOV
** Objective Magnification

