

SAFETY DATA SHEET

Version 1

Section 1: Identification of Product and Company

1.1 Product Identifiers

Product Name: **CLARITY Lipid Clearing Solution**
Manufacturer: ClearLight Biotechnologies

1.2 Relevant identified uses of substance and uses advised against

Identified uses: Laboratory reagent, For Research Use Only.

1.3 Details of the supplier of the safety data sheet

Company: ClearLight Biotechnologies, Inc.
428 Oakmead Parkway
Sunnyvale, CA 94085
Telephone: (800) 251-8905
Email address: contact@clearlightbiotech.com

1.4 Emergency telephone number

Emergency phone: (800) 251-8905

Section 2: Hazards identification

2.1 Classification of substance

GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Reproductive Toxicity (Category 1B), H360

2.2 GHS Label elements, including precautionary statements



Pictogram

Signal Word: **Danger**

Hazard Statement(s):

H360: May damage fertility or the unborn child.

Prevention Statement(s):

P201: Obtain special instructions before use.



428 Oakmead Parkway
Sunnyvale, CA 94085
Main: (800) 251-8905
URL: www.clearlightbiotechnologies.com
Email: laborders@clearlightbiotech.com

P202: Do not handle until all safety precautions have been read and understood.

P281: Use Personal protective equipment as required.

Response Statement(s):

P308+P313: IF exposed or concerned: Get medical advice/attention.

Storage Statement(s):

P405: Store locked up.

Disposal Statement(s):

P501: Dispose of contents/container in accordance with local/regional/national/international regulations.

2.3 Hazards not otherwise classified (HNOC) or not covered by GHS – none

Section 3: Composition/information on ingredients

3.1 Substance

Component	CAS Number	Concentration
1-Thioglycerol	96-27-5	0-40% ¹
20X Borate Buffer	10043-35-3	0-40% ²
20% SDS	151-21-3	0-40% ³

^{1,2} - Exact concentrations withheld as trade secrets

Section 4: First Aid measures

4.1 Description of first aid measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

Consult a physician.

In case of skin contact

Wash off with soap and plenty of water. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Consult a physician.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available.

Section 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

5.2 Special hazards arising from the substance or mixture

Carbon oxides, Sulphur oxides

5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

5.4 Further Information

No data available.

Section 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Use personal protective equipment. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas.

For personal protection see section 8.

6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

Discharge into the environment must be avoided.

6.3 Methods and materials for containment and cleaning up

Soak up with inert absorbent material and dispose of as hazardous waste. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

For disposal see section 13.



428 Oakmead Parkway
Sunnyvale, CA 94085
Main: (800) 251-8905
URL: www.clearlightbiotechnologies.com
Email: laborers@clearlightbiotech.com

Section 7: Handling and storage

7.1 Precautions for safe handling

Avoid contact with skin and eyes. Avoid inhalation of vapor or mist. Wash skin after handling. Do not eat, drink or smoke when using this product. Avoid release to the environment. Wear protective gloves. Take off contaminated clothing and wash before reuse. Collect spillage.

7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

7.3 Specific end use(s)

Aside from the uses mentioned in section 1.2, no other specific uses are specified.

Section 8: Exposure controls/personal protection

8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Appropriate engineering controls

Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product.

Personal protective equipment

Skin protection

Handle with gloves and dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands after use or handling.

Recommended glove material: Nitrile rubber

Advisory recommendation only.

Eye protection

Safety goggles or glasses.

Body Protection

Lightweight protective clothing.

Control of environmental exposure

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

Discharge into the environment must be avoided.

Section 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance:	Clear liquid
Odor:	No data available
Odor Threshold:	No data available
pH as supplied:	8.5 at 25°C
Melting point:	No data available
Freezing point:	No data available
Boiling point:	No data available
Flash point:	No data available
Evaporation rate:	No data available
Flammability (solid, gas):	No data available
Upper/lower flammability or explosive limits:	No data available
Vapor pressure:	No data available
Vapor density:	No data available
Relative density:	No data available
Water solubility:	No data available
Partition coefficient:	No data available
Auto-ignition temperature:	No data available
Decomposition temperature:	No data available
Viscosity:	No data available
Explosive properties:	No data available
Oxidizing properties:	No data available

9.2 Other safety information

No data available.

Section 10: Stability and reactivity

10.1 Reactivity

No data available.

10.2 Chemical stability

Stable under recommended storage conditions.



428 Oakmead Parkway
Sunnyvale, CA 94085
Main: (800) 251-8905
URL: www.clearlightbiotechnologies.com
Email: laborers@clearlightbiotech.com

10.3 Possibility of hazardous reactions

No data available.

10.4 Conditions to avoid

No data available.

10.5 Incompatible materials

Strong bases, Strong oxidizing agents.

10.6 Hazardous decomposition products

Hazardous decomposition products formed under fire conditions: Carbon oxides.

Other decomposition products: No data available.

In the event of fire: see section 5.

Section 11: Toxicological information

11.1 Information on toxicological effects

Acute toxicity

Product Information

Oral LD50

Based on ATE data, the classification criteria are not met. ATE > 5,000 mg/kg.

Dermal LD50

Based on ATE data, the classification criteria are not met. ATE > 5,000 mg/kg.

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitization

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

IARC:

No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

NTP:

No component of this product present at levels greater than or



428 Oakmead Parkway
Sunnyvale, CA 94085
Main: (800) 251-8905
URL: www.clearlightbiotechnologies.com
Email: laborers@clearlightbiotech.com

equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

Reproductive toxicity

No data available

Specific target organ toxicity - single exposure

No data available

Specific target organ toxicity - repeated exposure

No data available

Aspiration hazard

No data available

Additional Information

RTECS: Not available

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12: Ecological information

12.1 Toxicity

No data available.

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not Conducted.

12.6 Other adverse effects



428 Oakmead Parkway
Sunnyvale, CA 94085
Main: (800) 251-8905
URL: www.clearlightbiotechnologies.com
Email: laborders@clearlightbiotech.com

Harmful to aquatic life. An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

Section 13: Disposal considerations

13.1 Waste treatment methods

Product

Offer surplus and non-recyclable solutions to a licensed disposal company. Contact a licensed professional waste disposal service to dispose of this material.

Section 14: Transport information

DOT (US)

Not dangerous goods.

IMDG

Not dangerous goods.

IATA

Not dangerous goods.

Section 15: Regulatory information

SARA 302 Components

This material does not contain any components with a section 302 EHS TPQ.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

Acute Health Hazard

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16: Other information

Further Information:

©2020 ClearLight Biotechnologies, Inc. License granted to make unlimited paper copies for internal use only. The above information is believed to be correct, but does not purport to be



428 Oakmead Parkway
Sunnyvale, CA 94085
Main: (800) 251-8905
URL: www.clearlightbiotechnologies.com
Email: laborders@clearlightbiotech.com

all inclusive and shall be used only as a guide. It is the responsibility of the user to determine the application of this data to the necessary safety precautions. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. ClearLight Biotechnologies, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. See www.clearlightbiotechnologies.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

The branding on the header of this document may temporarily not visually match the product purchased in cases of brand transition. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact us at contact@clearlightbiotech.com.

Date of Preparation/Last Revision: October 7, 2020
